

E 451

Mesh trays

for optimum accommodation of different instruments.



- Mesh tray with lid for small pieces
- Wire mesh: Base 1 mm, sides 0.8 mm, lid 1 mm
- Mesh width: Base 3 mm, sides 1.7 mm, lid 3 mm

EAN: 4002511908154 / Material number: 04203910 / Old Material Number: 69545101D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Mesh tray	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of anesthetic instruments	•
Reprocessing of MIS instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	2 3/16 (55) in (mm)
External dimensions, net width in inches	6 5/16 (160) in (mm)
External dimensions, net depth in inches	9 1/2 (240) in (mm)
External dimensions, gross height in inches	3 (75)
External dimensions, gross width in inches	6 1/2 (165)
External dimensions, gross depth in inches	9 11/16 (245)
Net weight in lbs	2 (0.79)
Gross weight in lbs	2 (0.85)