

E 147/1

Mesh Basket

for the optimum loading of 10-12 mouthwash glasses.



- For 10-12 mouth rinse tumblers, max. Ø 80 mm
- Plastic-coated

EAN: 4002512122313 / Material number: 04746620 / Old Material Number: 68514703D

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
Product category	
Other inserts	•
Product characteristics	
Material	Wire, rilsan-coated
Color	White
Application area	
Reprocessing of dental instruments	•
Capacity	
Mouth rinse cups [number]	12
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	155
External dimensions, net width in inches	220
External dimensions, net depth in inches	455
External dimensions, gross height in inches	10 1/4 (260)
External dimensions, gross width in inches	20 1/8 (510)
External dimensions, gross depth in inches	9 1/2 (240)
Net weight in lbs	3 (1.29)
Gross weight in lbs	4 (1.59)