

E 468

Mesh Basket

for the optimum loading of different utensils.



- Made of welded stainless steel mesh with lid and sub-divider
- Mesh width 5 x 5 x 1 mm

EAN: 4002512277280 / Material number: 05043640 / Old Material Number: 69546801D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Mesh tray	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of anesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	76
External dimensions, net width in inches	170
External dimensions, net depth in inches	250
External dimensions, gross height in inches	4 (100)
External dimensions, gross width in inches	8 1/8 (205)
External dimensions, gross depth in inches	12 1/2 (315)
Net weight in lbs	2 (0.67)
Gross weight in lbs	2 (0.82)