

AUF 1 Holder

for optimum reprocessing of transmission instruments in upper basket O 177/1.



- Comprises: seal, insert nut, clamp
- Rack lower part, rack upper part, 5 filter plates

EAN: 4002512427517 / Material number: 05159510 / Old Material Number: 68551101D

Product category	
Thermal disinfectors	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8591	•
Product category	
Holder and/or adapter for transmission instruments	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of dental instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	0 (0) in (mm)
External dimensions, net width in inches	0 (0) in (mm)
External dimensions, net depth in inches	0 (0) in (mm)
External dimensions, gross height in inches	3 3/16 (80)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	3 15/16 (100)
Net weight in lbs	2 (0.90)
Gross weight in lbs	5 (1.84)
Accessory pack	
FP filter plates	•