

ML
Magnets
for automatic recognition of mobile unit/spray arm monitoring



- 15 different configurable codes
- 6 magnets for the magnetic strip installed for each area of use

EAN: 4002512423083 / Material number: 05267440 / Old Material Number: 69745702D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
Product category	
Magnetic strip for automatic mobile unit recognition	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel, Black
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	0
External dimensions, net width in inches	0
External dimensions, net depth in inches	0
External dimensions, gross height in inches	3 (75)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	2 13/16 (70)
Net weight in lbs	2 (0.85)
Gross weight in lbs	3 (1.00)