

ML

Magnets

for automatic recognition of mobile unit/spray arm monitoring



- 15 different configurable codes
- 6 magnets for the magnetic strip installed for each area of use

EAN: 4002512423083 / Material number: 05267440 / Old Material Number: 69745702D

<b>Product category</b>	
Large-chamber laboratory glasswashers	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
<b>Product category</b>	
Magnetic strip for automatic mobile unit recognition	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel, Black
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, gross height in inches	3 (75)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	2 13/16 (70)
Net weight in lbs	2 (0.85)
Gross weight in lbs	2 (1.00)