

E 520
 Mesh Basket
 for the optimum loading of root canal instruments.



- Suitable for 17 root canal instruments
- Hinged
- Can be steam sterilised at 121 °C/134 °C

EAN: 4002512472524 / Material number: 05409080 / Old Material Number: 68552001D

Product category	
Thermal disinfectors	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
Product category	
Other inserts	•
Product characteristics	
Material	Plastic
Color	White
Capacity	
Root-canal instruments [number]	17
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	45
External dimensions, net width in inches	75
External dimensions, net depth in inches	30
External dimensions, gross height in inches	3 1/4 (80)
External dimensions, gross width in inches	3 3/4 (95)
External dimensions, gross depth in inches	4 3/8 (110)
Net weight in lbs	1 (0.05)
Gross weight in lbs	1 (0.10)