


E 747

Module

for optimum loading of centrifuge tubes, phials and test tubes, for example.



- 104 nozzles for centrifuge tubes, phials, test tubes, etc.
- 104 nozzles 2.5 x 110 mm

EAN: 4002512591218 / Material number: 05464630 / Old Material Number: 69574701D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
<b>Product category</b>	
Insert/module for test tubes	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of pipettes/vials	•
<b>Capacity</b>	
Test tubes [number]	104
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	6 5/8 (168)
External dimensions, net width in inches	19 3/8 (492)
External dimensions, net depth in inches	19 9/16 (496)
External dimensions, gross height in inches	25 5/8 (650)
External dimensions, gross width in inches	26 (660)
External dimensions, gross depth in inches	31 15/16 (810)
Net weight in lbs	84 (38.00)
Gross weight in lbs	98 (44.00)