


E 523

Mesh Basket

for the optimum loading of 6 or more cassettes.



- Equipped with 7 holders (6 compartments), H 145, spacing 50 mm

EAN: 4002513199116 / Material number: 05649700 / Old Material Number: 69552301D

|   |                 |
|---|-----------------|
| <b>Product category</b>                     |                 |
| Thermal disinfectors                        | •               |
| <b>Machine affiliation</b>                  |                 |
| PG 8581                                     | •               |
| PG 8591                                     | •               |
| <b>Product category</b>                     |                 |
| Insert for trays                            | •               |
| <b>Product characteristics</b>              |                 |
| Material                                    | Stainless steel |
| Color                                       | Stainless Steel |
| Number of compartments                      | 6               |
| <b>Application area</b>                     |                 |
| Reprocessing of dental instruments          | •               |
| <b>Standard electrical connection</b>       |                 |
| Number of phases                            | 0               |
| Voltage in V                                | 0               |
| Frequency in Hz                             | 0               |
| Total rated load in kW                      | 0.00            |
| Fuse rating in A                            | 0               |
| <b>Dimensions and weight</b>                |                 |
| External dimensions, net height in inches   | 151             |
| External dimensions, net width in inches    | 220             |
| External dimensions, net depth in inches    | 450             |
| External dimensions, gross height in inches | 7 3/4 (195)     |
| External dimensions, gross width in inches  | 9 1/2 (240)     |
| External dimensions, gross depth in inches  | 19 1/2 (495)    |
| Net weight in lbs                           | 4 (1.40)        |
| Gross weight in lbs                         | 4 (1.65)        |