

E 476

Holder

for micro instruments with a Ø of 4-8 mm.



- For use in trays with 5 mm mesh
- 50 items per bag

EAN: 4002513239454 / Material number: 05692400 / Old Material Number: 69747601D

Product category

Thermal disinfectors	•
Cleaning machine, small, dental	•
SlimLine washer-disinfectors, dental	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•

Machine affiliation

G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWD 8531	•
PWD 8532	•
PWA 8672	•
PWD 7121	•
PWD 7122	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•

Product category

Holder	•
--------	---

Application area

Reprocessing of MIS instruments	•
Reprocessing of ophthalmic instruments	•
Reprocessing of dental instruments	•

Product characteristics

Material	Plastic
Color	Gray

Standard electrical connection

Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0

Dimensions and weight

External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)

E 476

Holder

for micro instruments with a Ø of 4-8 mm.

EAN: 4002513239454 / Material number: 05692400 / Old Material Number: 69747601D



External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	0 (0.08)
Gross weight in lbs	0 (0.08)