


E 470

Injector nozzle

for use in the mobile injector unit, Ø 2.5 x 90 mm.

EAN: 4002513248814 / Material number: 05701580 / Old Material Number: 69547001D



- Suitable for mobile injector units, with screw connection
- Suitable for bottles, Erlenmeyer flasks and round bottomed flasks for example

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Wire, rilsan-coated, Stainless steel
Color	Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	3 9/16 (90)
External dimensions, net width in inches	0 1/8 (3)
External dimensions, net depth in inches	0 1/8 (3)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 15/16 (100)
External dimensions, gross depth in inches	5 15/16 (150)
Diameter in inches	0 1/8 (3)
Net weight in lbs	1 (0.02)
Gross weight in lbs	1 (0.03)