


E 755

Module

for the optimum loading of narrow neck- glassware 25 - 100 ml.



- Equipped with 36 nozzles (E 470) 2.5 x 90 mm with holder

EAN: 4002513248821 / Material number: 05701590 / Old Material Number: 69575501D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
G 7826	•
<b>Product category</b>	
Insert/module for narrow-neck glassware	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of pipettes/vials	•
<b>Capacity</b>	
Narrow-neck glassware [number]	36
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	5 1/8 (130)
External dimensions, net width in inches	19 3/8 (492)
External dimensions, net depth in inches	19 9/16 (496)
External dimensions, gross height in inches	21 1/2 (545)
External dimensions, gross width in inches	6 5/16 (160)
External dimensions, gross depth in inches	20 7/8 (530)
Net weight in lbs	11 (4.80)
Gross weight in lbs	12 (5.40)