

ML/2

Magnetic strip

for automatic recognition of mobile unit/spray arm monitoring



- ML/2 is required for spray arm monitoring
- 5 magnets, individually configurable
- 15 possible combinations

EAN: 4002513274158 / Material number: 05723270 / Old Material Number: 69745703D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Magnetic strip for automatic mobile unit recognition	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel, Black
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	0 3/8 (9) in (mm)
External dimensions, net width in inches	1 (25) in (mm)
External dimensions, net depth in inches	10 7/8 (275) in (mm)
External dimensions, gross height in inches	3 3/4 (95)
External dimensions, gross width in inches	8 1/8 (205)
External dimensions, gross depth in inches	12 7/16 (315)
Net weight in lbs	1 (0.17)
Gross weight in lbs	1 (0.32)