

E 957

Mobile unit

for the optimum loading of up to 12 large volume laboratory glassware.



- Equipped with 12 nozzles with Ø 10 mm
- Height adjustable holding frame
- Magnetic strip for automatic mobile unit recognition
- Max. vertical clearance from star overlay: 24 1/4 (615) in (mm) mm

EAN: 4002513298789 / Material number: 05746300 / Old Material Number: 69595701D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8527	•
<b>Product category</b>	
Other mobile units	•
Module for large-volume laboratory glassware	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Maximum vertical clearance in mm	24 1/4 (615)
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Large-volume laboratory glassware [number]	12
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	13 15/16 (353)
External dimensions, net width in inches	25 1/4 (640)
External dimensions, net depth in inches	31 1/8 (790)
External dimensions, gross height in inches	15 3/8 (390)
External dimensions, gross width in inches	26 3/16 (665)
External dimensions, gross depth in inches	31 3/4 (805)
Net weight in lbs	38 (16.98)
Gross weight in lbs	39 (17.45)