

E 131/1

Mesh Basket

for optimum loading of 5 cassettes or kidney trays.



- Suitable for 5 mesh or kidney trays
- Equipped with 6 holders, height: 6 5/16 (160) in (mm) mm, distance: 3 3/16 (80) in (mm) mm

EAN: 4002513317336 / Material number: 05771030 / Old Material Number: 69513102D

| | |
|---|-----------------|
| Product category | |
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8535 | • |
| PG 8536 | • |
| PG 8581 | • |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8591 | • |
| PG 8592 | • |
| Product category | |
| Insert for mesh trays | • |
| Insert for kidney bowls/basins | • |
| Product characteristics | |
| Material | Stainless steel |
| Color | Stainless Steel |
| Distance between rails in mm | 3 3/16 (80) |
| Rail height in mm | 6 5/16 (160) |
| Application area | |
| Reprocessing of kidney bowls | • |
| Reprocessing of surgical instruments | • |
| Capacity | |
| Kidney bowls [number] | 5 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0-0 |
| Frequency in Hz | 0-0 |
| Total rated load in kW | 0.00-0.00 |
| Fuse rating in A | 0-0 |
| Dimensions and weight | |
| External dimensions, net height in inches | 6 5/8 (168) |
| External dimensions, net width in inches | 7 1/8 (180) |
| External dimensions, net depth in inches | 19 1/2 (495) |
| External dimensions, gross height in inches | 7 5/16 (185) |
| External dimensions, gross width in inches | 7 7/8 (200) |
| External dimensions, gross depth in inches | 20 7/8 (530) |
| Net weight in lbs | 2 (0.61) |
| Gross weight in lbs | 2 (0.90) |