

E 131/1

Mesh Basket

for optimum loading of 5 cassettes or kidney trays.



- Suitable for 5 mesh or kidney trays
- Equipped with 6 holders, height: 6 3/8 (160) in (mm) mm, distance: 3 1/4 (80) in (mm) mm

EAN: 4002513317336 / Material number: 05771030 / Old Material Number: 69513102D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Insert for mesh trays	•
Insert for kidney bowls/basins	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Distance between rails in mm	3 1/4 (80)
Rail height in mm	6 3/8 (160)
Application area	
Reprocessing of kidney bowls	•
Reprocessing of surgical instruments	•
Capacity	
Kidney bowls [number]	5
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	168
External dimensions, net width in inches	180
External dimensions, net depth in inches	495
External dimensions, gross height in inches	7 3/8 (185)
External dimensions, gross width in inches	7 7/8 (200)
External dimensions, gross depth in inches	20 7/8 (530)
Net weight in lbs	2 (0.61)
Gross weight in lbs	2 (0.90)