


E 965  
 Mesh Basket  
 for the optimum loading of beakers 250 - 600 ml.



- Equipped with 15 x 3 stay bars for beakers from 250 ml to 600 ml

EAN: 4002513369229 / Material number: 05848290 / Old Material Number: 69596501D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8527	•
<b>Product category</b>	
Insert/module for beakers	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Beakers [number]	98
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	173
External dimensions, net width in inches	357
External dimensions, net depth in inches	522
External dimensions, gross height in inches	11 1/16 (280)
External dimensions, gross width in inches	21 5/16 (540)
External dimensions, gross depth in inches	22 1/16 (560)
Net weight in lbs	5 (2.07)
Gross weight in lbs	7 (2.86)