


E 963

Mesh Basket

for the optimum loading of beakers up to 250 ml.



- Equipped with 33 x 3 stay bars for beakers up to 250 ml

EAN: 4002513369236 / Material number: 05848300 / Old Material Number: 69596301D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8527	•
<b>Product category</b>	
Insert/module for beakers	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Beakers [number]	98
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	155
External dimensions, net width in inches	357
External dimensions, net depth in inches	522
External dimensions, gross height in inches	11 1/4 (285)
External dimensions, gross width in inches	21 3/4 (550)
External dimensions, gross depth in inches	20 7/8 (530)
Net weight in lbs	4 (1.58)
Gross weight in lbs	4 (1.70)