


E 960/1

Mesh Basket

for the optimum loading of wide-neck Erlenmeyer flasks.



- Equipped with 20 large and 26 small spring hooks
- Suitable for Erlenmeyer flasks & measuring cylinders

EAN: 4002513400953 / Material number: 05892360 / Old Material Number: 69596002D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
Product category	
Insert/module for wide-neck glassware	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Erlenmeyer flasks [number]	46
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	185
External dimensions, net width in inches	357
External dimensions, net depth in inches	522
External dimensions, gross height in inches	11 1/16 (280)
External dimensions, gross width in inches	20 7/8 (530)
External dimensions, gross depth in inches	22 1/16 (560)
Net weight in lbs	2 (0.80)
Gross weight in lbs	2 (0.90)