

E 441/1

Mesh Basket

for the optimum loading of Micro instruments.



- The mesh width on the base is 1.7 mm
- The sides are closed, stackable
- Sections w/ 6 adjustable dividers for positioning of instruments
- Can be steam sterilised at 121 °C/134 °C

EAN: 4002513460476 / Material number: 05988240 /
Old Material Number: 69544102D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
SlimLine washer-disinfectors, dental	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 7121	•
PWD 7122	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert/module for MIS instruments	•
Application area	
Reprocessing of MIS instruments	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	2 3/8 (60)
External dimensions, net width in inches	7 1/4 (183)
External dimensions, net depth in inches	11 3/16 (284)
External dimensions, gross height in inches	3 15/16 (100)
External dimensions, gross width in inches	8 5/16 (210)
External dimensions, gross depth in inches	12 7/16 (315)
Net weight in lbs	2 (0.77)
Gross weight in lbs	2 (0.92)