

E 441/1

Mesh Basket

for the optimum loading of Micro instruments.



- The mesh width on the base is 1.7 mm
- The sides are closed, stackable
- Sections w/ 6 adjustable dividers for positioning of instruments
- Can be steam sterilised at 121 °C/134 °C

EAN: 4002513460476 / Material number: 05988240 / Old Material Number: 69544102D

| | |
|---|--------------------------|
| Product category | |
| Thermal disinfectors | • |
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| G 7831 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8581 | • |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8591 | • |
| PG 8592 | • |
| Product category | |
| Insert/module for MIS instruments | • |
| Product characteristics | |
| Material | Stainless steel, Plastic |
| Color | Stainless Steel, Black |
| Application area | |
| Reprocessing of MIS instruments | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0-0 |
| Frequency in Hz | 0-0 |
| Total rated load in kW | 0.00-0.00 |
| Fuse rating in A | 0-0 |
| Dimensions and weight | |
| External dimensions, net height in inches | 2 3/8 (60) |
| External dimensions, net width in inches | 7 1/4 (183) |
| External dimensions, net depth in inches | 11 3/16 (284) |
| External dimensions, gross height in inches | 3 15/16 (100) |
| External dimensions, gross width in inches | 8 5/16 (210) |
| External dimensions, gross depth in inches | 12 7/16 (315) |
| Net weight in lbs | 2 (0.77) |
| Gross weight in lbs | 3 (0.92) |