


U 184/1

Lower basket

for the optimum loading of centrifuge tubes, phials and test tubes.



- Equipped with 96 injector nozzles
- Only suitable for use with upper baskets O 175 / O 187 / O 184
- Vertical clearance: Approx. 6 3/4 (170) in (mm) mm

EAN: 4002513805659 / Material number: 06439180 / Old Material Number: 69118403D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8536	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Maximum vertical clearance in mm	6 3/4 (170)
Application area	
Reprocessing of pipettes/vials	•
Reprocessing of test tubes	•
Capacity	
Test tubes [number]	96
Phials [number]	96
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	9 7/8 (250)
External dimensions, net width in inches	20 5/16 (515)
External dimensions, net depth in inches	20 15/16 (531)
External dimensions, gross height in inches	10 11/16 (270)
External dimensions, gross width in inches	20 11/16 (525)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	12 (5.00)
Gross weight in lbs	14 (6.00)
Accessory pack	
ID 90, injector nozzle, Ø 1/8" (2.5 mm), length 3 9/16" (90 mm) [number]	96