


U 874/1

Lower basket

for the optimum loading of different inserts and mesh trays.



- Vertical clearance can vary with different upper baskets
- With bracket for magnetic strip for ML/2 autom. unit recognition

EAN: 4002513808919 / Material number: 06442940 / Old Material Number: 69187403D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Carrier/basket carrier	•
Lower basket	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of MIS instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of medical shoes	•
Reprocessing of basins	•
Reprocessing of kidney bowls	•
Reprocessing of surgical instruments	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	2 (50)
External dimensions, net width in inches	20 15/16 (531)
External dimensions, net depth in inches	20 5/16 (515)
External dimensions, gross height in inches	3 15/16 (100)
External dimensions, gross width in inches	20 7/8 (530)
External dimensions, gross depth in inches	21 1/2 (545)
Net weight in lbs	5 (2.10)
Gross weight in lbs	7 (2.97)