


A 12/1

Mesh Basket

For the gentle cleaning of delicate items by interrupting the spray jet.



- With perforations 7 x 7 mm
- For use in upper and lower basket
- Recommended for use in PG 8536

EAN: 4002513975697 / Material number: 06629670 / Old Material Number: 69501203D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Spacer	•
Retaining frame	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	0 1/16 (1) in (mm)
External dimensions, net width in inches	8 7/8 (224) in (mm)
External dimensions, net depth in inches	16 15/16 (429) in (mm)
External dimensions, gross height in inches	3 (75)
External dimensions, gross width in inches	11 1/16 (280)
External dimensions, gross depth in inches	20 1/8 (510)
Net weight in lbs	1 (0.37)
Gross weight in lbs	2 (0.56)