

E 430/1

Mesh Basket

for optimum accommodation of different instruments.



- Comprising wire mesh with 5 mm mesh width

EAN: 4002513984422 / Material number: 06643200 / Old Material Number: 69743003D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Mesh tray	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	40
External dimensions, net width in inches	150
External dimensions, net depth in inches	445
External dimensions, gross height in inches	3 15/16 (100)
External dimensions, gross width in inches	6 5/16 (160)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	1 (0.36)
Gross weight in lbs	2 (0.48)