


E 522/1

Mesh Basket

for the optimum loading of 9 impression trays.



- Equipped with 9 holders

EAN: 4002514072722 / Material number: 06697010 / Old Material Number: 68552202D

Product category	
Thermal disinfectors	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
Product category	
Insert/module for dental instruments	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of dental instruments	•
Capacity	
Impression tray [number]	9
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	140
External dimensions, net width in inches	100
External dimensions, net depth in inches	190
External dimensions, gross height in inches	6 3/4 (170)
External dimensions, gross width in inches	9 3/4 (245)
External dimensions, gross depth in inches	9 1/2 (240)
Net weight in lbs	1 (0.26)
Gross weight in lbs	1 (0.36)