

E 404/1  
Mobile unit  
for the optimum loading of up to 38 pipettes.



- Suitable for 38 pipettes in 3 rows
- 1st Row: 10 pipettes x 100 ml (up to 550 mm), bar distance 20 mm
- 2nd Row: 14 pipettes x 25 ml, bar distance 26 mm
- 3rd Row: 14 pipettes x 10 ml, bar distance 26 mm
- Without bracket for magnetic strip for mobile unit recognition

EAN: 4002514087061 / Material number: 06727910 / Old Material Number: 69540402D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Mobile unit for pipettes	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Number of rows	3
<b>Application area</b>	
Reprocessing of pipettes/vials	•
<b>Capacity</b>	
Pipettes [number]	38
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	19 13/16 (502)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	20 3/8 (517)
External dimensions, gross height in inches	19 11/16 (500)
External dimensions, gross width in inches	21 11/16 (550)
External dimensions, gross depth in inches	21 11/16 (550)
Net weight in lbs	13 (5.70)
Gross weight in lbs	17 (7.40)