

E 405/1

Mobile unit

for optimum loading of up to 38 pipettes, connection for hot air drying.



- Suitable for 38 pipettes in 3 rows
- 1st Row: 10 pipettes x 100 ml (up to 550 mm), bar distance 20 mm
- 2nd Row: 14 pipettes x 25 ml, bar distance 26 mm
- 3rd Row: 14 pipettes x 10 ml, bar distance 26 mm
- Bracket for ML/2 for mobile unit recognition / spray arm control

EAN: 4002514087078 / Material number: 06727920 / Old Material Number: 69540502D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Mobile unit with connection for hot-air drying	•
Mobile unit for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Number of rows	3
Application area	
Reprocessing of pipettes/vials	•
Capacity	
Pipettes [number]	38
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	19 13/16 (502)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	20 3/8 (517)
External dimensions, gross height in inches	19 11/16 (500)
External dimensions, gross width in inches	21 5/16 (540)
External dimensions, gross depth in inches	20 7/8 (530)
Net weight in lbs	14 (6.20)
Gross weight in lbs	18 (8.10)