


E 105/1

Mesh Basket

for the optimum loading of test tubes up to 12 x 165 mm.

EAN: 4002514231075 / Material number: 06907650 / Old Material Number: 69510502D



- Suitable for about 200 Test tubes (up to 12 x 165 mm)
- Capacity including lid A 13
- The mesh width is 9 x 9 mm

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Insert/module for test tubes	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Number of compartments	6
Application area	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	200
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	8 3/8 (212)
External dimensions, net width in inches	7 7/8 (200)
External dimensions, net depth in inches	8 5/16 (210)
External dimensions, gross height in inches	9 5/16 (235)
External dimensions, gross width in inches	9 1/2 (240)
External dimensions, gross depth in inches	9 7/8 (250)
Net weight in lbs	2 (0.90)
Gross weight in lbs	3 (1.14)
Accessory pack	
A 13 lid [number]	1