

E 139/1

Mesh Basket

for the optimum loading of test tubes up to 12 x 200 mm.



- Suitable for test tubes up to 12 x 200 mm
- Divided into 6 sections
- Mesh width 9 x 9 mm
- For use in lower basket
- Including lid A 13

EAN: 4002514231082 / Material number: 06907660 / Old Material Number: 69513902D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Insert/module for test tubes	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Number of compartments	6
Application area	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	200
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	9 5/8 (243)
External dimensions, net width in inches	7 7/8 (200)
External dimensions, net depth in inches	9 1/16 (230)
External dimensions, gross height in inches	9 11/16 (245)
External dimensions, gross width in inches	9 11/16 (245)
External dimensions, gross depth in inches	10 7/16 (265)
Net weight in lbs	3 (1.26)
Gross weight in lbs	3 (1.29)
Accessory pack	
A 13 lid [number]	1