


O 188/2

Upper basket

for the optimum loading of different inserts.



- Vertical clearance: 165 mm, height adjustable +/- 20 mm
- With integrated magnetic spray arm
- Magnet for use with spray arm sensors for PG 8535, PG 8536

EAN: 4002514556178 / Material number: 07333030 / Old Material Number: 69118802D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Carrier/basket carrier	•
Upper basket	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	11 1/16 (280)
External dimensions, net width in inches	8 1/2 (215)
External dimensions, net depth in inches	18 3/4 (475)
External dimensions, gross height in inches	11 1/4 (285)
External dimensions, gross width in inches	21 5/16 (540)
External dimensions, gross depth in inches	21 7/8 (555)
Net weight in lbs	6 (2.40)
Gross weight in lbs	7 (3.16)