

O 190/2
Upper basket
for the optimum loading of different inserts.



- Vertical clearance: 8 1/2 (215) in (mm) mm
- With integrated magnetic spray arm
- Magnet for use with spray arm sensors for PG 8535, PG 8536

EAN: 4002514556185 / Material number: 07333050 / Old Material Number: 69119002D

| | |
|---|-----------------|
| Product category | |
| Washer-disinfectors, laboratory | • |
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8535 | • |
| PG 8536 | • |
| Product category | |
| Carrier/basket carrier | • |
| Upper basket | • |
| Product characteristics | |
| Material | Stainless steel |
| Color | Stainless Steel |
| Maximum vertical clearance in mm | 8 1/2 (215) |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0-0 |
| Frequency in Hz | 0-0 |
| Total rated load in kW | 0.00-0.00 |
| Fuse rating in A | 0-0 |
| Dimensions and weight | |
| External dimensions, net height in inches | 10 7/16 (265) |
| External dimensions, net width in inches | 20 15/16 (531) |
| External dimensions, net depth in inches | 18 3/4 (475) |
| External dimensions, gross height in inches | 11 1/4 (285) |
| External dimensions, gross width in inches | 21 5/16 (540) |
| External dimensions, gross depth in inches | 21 7/8 (555) |
| Net weight in lbs | 6 (2.40) |
| Gross weight in lbs | 8 (3.19) |