

E 942/3

Module

for the optimum loading of up to 116 pipettes.



- Suitable for 116 pipettes up to 580 mm in height
- Compartment size ?? x ??

EAN: 4002514605036 / Material number: 07459390 / Old Material Number: 69594204D

|  |                 |
|--|-----------------|
| <b>Product category</b>                        |                 |
| Large-capacity washer-disinfectors, laboratory | •               |
| <b>Machine affiliation</b>                     |                 |
| PG 8527  | •               |
| PG 8528  | •               |
| <b>Product category</b>                        |                 |
| Insert/module for pipettes                     | •               |
| <b>Product characteristics</b>                 |                 |
| Material                                       | Stainless steel |
| Color  | Stainless Steel |
| <b>Application area</b>                        |                 |
| Reprocessing of pipettes/vials                 | •               |
| <b>Capacity</b>                                |                 |
| Pipettes [number]                              | 116             |
| <b>Standard electrical connection</b>          |                 |
| Number of phases                               | 0               |
| Voltage in V                                   | 0               |
| Frequency in Hz                                | 0               |
| Total rated load in kW                         | 0.00            |
| Fuse rating in A                               | 0               |
| <b>Dimensions and weight</b>                   |                 |
| External dimensions, net height in inches      | 11 (279)        |
| External dimensions, net width in inches       | 22 (558)        |
| External dimensions, net depth in inches       | 13 7/8 (352)    |
| External dimensions, gross height in inches    | 16 9/16 (420)   |
| External dimensions, gross width in inches     | 21 11/16 (550)  |
| External dimensions, gross depth in inches     | 22 1/2 (570)    |
| Net weight in lbs                              | 78 (35.00)      |
| Gross weight in lbs                            | 100 (45.00)     |