

E 944/2

Module

for the optimum loading of narrow-neck glassware 500 - 1,000 ml.



- 15 nozzles (E 352) 6 x 220 mm with clips (E 354)

EAN: 4002514605050 / Material number: 07459410 / Old Material Number: 69594403D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Insert/module for narrow-neck glassware	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	15
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	9 7/8 (250)
External dimensions, net width in inches	22 (558)
External dimensions, net depth in inches	13 7/8 (352)
External dimensions, gross height in inches	7 1/8 (180)
External dimensions, gross width in inches	21 11/16 (550)
External dimensions, gross depth in inches	22 1/2 (570)
Net weight in lbs	9 (3.92)
Gross weight in lbs	10 (4.42)