

E 521/2

Mesh Basket

for the optimum loading of 7 extraction, orthodontic and/or technical forceps.



- Compartment size: 21 x 80 mm
- Also suitable for forceps, such as those used in podiatry

EAN: 4002514616438 / Material number: 07476430 / Old Material Number: 68552103D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Insert/module for dental instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Capacity</b>	
Forceps [number]	7
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	5 3/8 (135)
External dimensions, net width in inches	3 15/16 (100)
External dimensions, net depth in inches	7 1/2 (189)
External dimensions, gross height in inches	4 3/4 (120)
External dimensions, gross width in inches	6 15/16 (175)
External dimensions, gross depth in inches	8 11/16 (220)
Net weight in lbs	1 (0.42)
Gross weight in lbs	2 (0.46)