

E 521/2

Mesh Basket

for the optimum loading of 7 extraction, orthodontic and/or technical forceps.

EAN: 4002514616438 / Material number: 07476430 / Old Material Number: 68552103D



- Compartment size: 21 x 80 mm
- Also suitable for forceps, such as those used in podiatry

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
Product category	
Insert/module for dental instruments	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of dental instruments	•
Capacity	
Forceps [number]	7
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	5 3/8 (135) in (mm)
External dimensions, net width in inches	3 15/16 (100) in (mm)
External dimensions, net depth in inches	7 1/2 (189) in (mm)
External dimensions, gross height in inches	4 3/4 (120)
External dimensions, gross width in inches	6 15/16 (175)
External dimensions, gross depth in inches	8 11/16 (220)
Net weight in lbs	1 (0.42)
Gross weight in lbs	2 (0.46)