

DOS S 20
 Bellows dispenser pump
 For dispensing liquid acid detergent.

EAN: 4002514641904 / Material number: 07514600 / Old Material Number: 69226004D



- For dispensing rinse aid, neutralizing agent
- For retro-fitting by Miele Service

| | |
|---|--------------|
| Product category | |
| Washer-disinfectors, laboratory | • |
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8536 | • |
| Product category | |
| Dispenser system | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0.00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in inches | 0 |
| External dimensions, net width in inches | 0 |
| External dimensions, net depth in inches | 0 |
| External dimensions, gross height in inches | 9 3/4 (245) |
| External dimensions, gross width in inches | 10 1/4 (258) |
| External dimensions, gross depth in inches | 10 3/4 (270) |
| Net weight in lbs | 6 (2.45) |
| Gross weight in lbs | 7 (2.91) |