

DOS NA 120

Bellows dispenser pump

For dispensing liquid neutral/alkaline detergent.



- For dispensing liquid cleaning agent, emulgator
- For dispensing chemical disinfectant
- For retro-fitting by Miele Service

EAN: 4002514641911 / Material number: 07514610 / Old Material Number: 69226005D

Product category	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8536	•
Product category	
Dispenser system	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	0
External dimensions, net width in inches	0
External dimensions, net depth in inches	0
External dimensions, gross height in inches	11 13/16 (300)
External dimensions, gross width in inches	9 7/8 (250)
External dimensions, gross depth in inches	16 9/16 (420)
Net weight in lbs	6 (2.45)
Gross weight in lbs	7 (2.89)