

E 491

Mesh Basket

for the optimum loading of rotating instruments.



- Suitable for 18 turbine shafts
- Suitable for 30 contra angle shafts
- Also suitable for rotating instruments, such as used in podiatry

EAN: 4002514697949 / Material number: 07600210 / Old Material Number: 68549101D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
Product category	
Insert/module for dental instruments	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of dental instruments	•
Capacity	
Sheaths [number]	48
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	2 1/8 (53)
External dimensions, net width in inches	2 13/16 (70)
External dimensions, net depth in inches	4 3/4 (120)
External dimensions, gross height in inches	3 15/16 (100)
External dimensions, gross width in inches	10 7/16 (265)
External dimensions, gross depth in inches	7 1/8 (180)
Net weight in lbs	2 (0.85)
Gross weight in lbs	3 (1.00)