

E 935/2

Mobile unit

for the optimum cleaning of instruments on 3 levels.



- Loading capacity: 9 DIN mesh trays\*
- Equipped with 2 integrated spray arms
- The loading areas for the individual levels are different
- Magnetic strip for mobile unit recognition (incl. magnets)

\* DIN mesh tray 485 x 254 x 50 mm

EAN: 4002514816432 / Material number: 07765780 / Old Material Number: 69593503D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
<b>Capacity</b>	
DIN mesh trays [number]	9
SPRI mesh trays [number]	6
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	20 11/16 (524)
External dimensions, net width in inches	25 1/4 (640)
External dimensions, net depth in inches	31 1/8 (790)
External dimensions, gross height in inches	28 3/8 (720)
External dimensions, gross width in inches	32 11/16 (830)
External dimensions, gross depth in inches	30 3/4 (780)
Net weight in lbs	59 (26.50)
Gross weight in lbs	67 (30.00)
<b>Accessory pack</b>	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1