


E 197

Mesh Basket

for optimum accommodation of different instruments.



- Mesh width on the base 3 mm, sides 1.7 mm, lid 8 mm
- Can be steam sterilised at 121 °C/134 °C

EAN: 4002514953380 / Material number: 07984090 / Old Material Number: 68519701D

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
<b>Product category</b>	
Mesh tray	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	42
External dimensions, net width in inches	150
External dimensions, net depth in inches	225
External dimensions, gross height in inches	3 (75)
External dimensions, gross width in inches	6 3/4 (170)
External dimensions, gross depth in inches	9 1/2 (240)
Net weight in lbs	2 (0.55)
Gross weight in lbs	2 (0.62)