


E 197

Mesh Basket

for optimum accommodation of different instruments.



- Mesh width on the base 3 mm, sides 1.7 mm, lid 8 mm
- Can be steam sterilised at 121 °C/134 °C

EAN: 4002514953380 / Material number: 07984090 / Old Material Number: 68519701D

Product category	
Thermal disinfectors	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mesh tray	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of dental instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	1 11/16 (42)
External dimensions, net width in inches	5 15/16 (150)
External dimensions, net depth in inches	8 7/8 (225)
External dimensions, gross height in inches	3 (75)
External dimensions, gross width in inches	6 3/4 (170)
External dimensions, gross depth in inches	9 1/2 (240)
Net weight in lbs	2 (0.55)
Gross weight in lbs	2 (0.62)