

E 337/1

Mesh Basket

for the optimum loading of upright instruments.



- Suitable for probes, tweezers or mouth specula
- Equipped with 4 plastic holders
- 12 compartments approx. width 22 x 28 mm
- 4 compartments approx. width 25 x 28 mm
- 48 compartments approx. width 13 x 14 mm

EAN: 4002514953397 / Material number: 07984130 / Old Material Number: 69533702D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
<b>Product category</b>	
Insert/module for dental instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	5 7/8 (148) in (mm)
External dimensions, net width in inches	6 13/16 (173) in (mm)
External dimensions, net depth in inches	17 9/16 (445) in (mm)
External dimensions, gross height in inches	8 1/8 (205)
External dimensions, gross width in inches	9 1/2 (240)
External dimensions, gross depth in inches	19 5/16 (490)
Net weight in lbs	5 (1.91)
Gross weight in lbs	5 (2.16)