

E 802/1

Mesh Basket

for the optimum loading of different upright instruments.



- Suitable for loading probes, tweezers, or mouth specula
- Equipped with 4 plastic holders
- 8 compartments approx. width 30 x 28 mm
- 8 compartments approx. width 16 x 28 mm
- 15 compartments approx. width 16 x 20 mm

EAN: 4002514953403 / Material number: 07984170 / Old Material Number: 69580202D

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
G 7831	•
PG 8581	•
PG 8562	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
<b>Product category</b>	
Insert for instruments	•
Insert/module for dental instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	148
External dimensions, net width in inches	163
External dimensions, net depth in inches	295
External dimensions, gross height in inches	7 1/2 (190)
External dimensions, gross width in inches	7 3/4 (195)
External dimensions, gross depth in inches	17 (430)
Net weight in lbs	2 (0.65)
Gross weight in lbs	3 (1.15)