


E 807

Mesh Basket

for optimum loading of 3 mesh or kidney trays.



- 4 holders (3 compartments), H 165 mm, spacing 2 11/16 (68) in (mm) mm

EAN: 4002514951638 / Material number: 07984220 / Old Material Number: 69580701D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
<b>Product category</b>	
Insert for mesh trays	•
Insert for kidney bowls/basins	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
Number of compartments	3
Distance between rails in mm	2 11/16 (68)
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Capacity</b>	
Diverse mesh trays [number]	3
Kidney bowls [number]	3
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	144
External dimensions, net width in inches	140
External dimensions, net depth in inches	290
External dimensions, gross height in inches	18 15/16 (480)
External dimensions, gross width in inches	12 5/8 (320)
External dimensions, gross depth in inches	5 1/8 (130)
Net weight in lbs	1 (0.31)
Gross weight in lbs	2 (0.47)