

DGG 20

Unperforated steam oven pan
for cooking food in gravy, stock, water (e.g. rice, pasta).

EAN: 4002515060650 / Material number: 08246340 / Old Material Number: 23996112D



- Capacity: 0.64 gallon, useable capacity: 0.48 gl
- Exterior dimensions W x D x H: 17.7 x 7.5 x 1.5 in
- Made from treated stainless steel, wall thickness .02 in
- Professional accessories manufactured in Europe
- Made to fit Miele Combi-Steam Ovens