

DGGL 12

Perforated steam oven pan

For blanching or cooking vegetables, fish, meat and potatoes and much more

EAN: 4002514849904 / Material number: 08249430 / Old Material Number: 23996108D



- Capacity: 1.43 gallon, useable capacity: 0.87 gal
- Exterior dimensions W x D x H: 17.6 x 15.2 x 1.57 in
- Made from treated stainless steel, wall thickness .02 in
- Professional accessories manufactured in Europe
- Made to fit Miele Combi-Steam Ovens