

DGGL 20

Perforated steam oven pan

For blanching or cooking vegetables, fish, meat and potatoes and much more

EAN: 4002515014509 / Material number: 08285410 / Old Material Number: 23996111D



- Capacity: 0.64 gallon, useable capacity: 0.48 gl
- Exterior dimensions W x D x H: 17.7 x 7.5 x 1.5 in
- Made from treated stainless steel, wall thickness .02 in
- Professional accessories manufactured in Europe
- Made to fit Miele Combi-Steam Ovens