

ProCare Universal 61  
 Regeneration salt, 3 x 2 kg (13.2 lb total)  
 for regeneration of the integrated water softener.



- Prevents calcium deposits and protects wash items

EAN: 4002515075319 / Material number: 09195780 / Old Material Number: 66916102EU1

Product category	
Fresh-water dishwashers	•
ProfiLine dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Product characteristics	
Medium type	Reactivation salt
Ingredients	Coarse grained salt
Prevents limescale and protects wash items	•
Dimensions and weight	
External dimensions, net height in inches	12 5/8 (320)
External dimensions, net width in inches	4 1/8 (104)
External dimensions, net depth in inches	8 15/16 (226)
External dimensions, gross height in inches	12 5/8 (320)
External dimensions, gross width in inches	3 15/16 (100)
External dimensions, gross depth in inches	8 7/8 (225)
Load size in g	211.64 (6000)
Gross weight in lbs	0.21 (6)
Compliance, test certificates, and symbols	
GHS/CLP compliant	•