

E 339/1

Mesh Basket

for the optimum loading of up to 12 tray dividers and trays.

EAN: 4002515094860 / Material number: 09233840 / Old Material Number: 69533902D



- Suitable for 12 tray dividers or trays
- Equipped with 13 holders (12 compartments).
- Maximum tray size 290 x 20 mm

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Insert for trays	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Number of compartments	12
Capacity	
Trays [number]	12
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	115
External dimensions, net width in inches	305
External dimensions, net depth in inches	498
External dimensions, gross height in inches	18 15/16 (480)
External dimensions, gross width in inches	12 5/8 (320)
External dimensions, gross depth in inches	5 1/8 (130)
Net weight in lbs	2 (0.80)
Gross weight in lbs	3 (1.22)